

Eliminate Bothersome Odors

Refresh Labo thoroughly eliminates annoying odors.

Also use Refresh Labo to wipe down tables and chairs and eliminate odors from furniture surfaces. Refresh Labo eliminates odors, germs, and viruses* and can even be used to wipe down surfaces, helping to simplify cleaning, reduce working times, and reduce the number of cleaning products you need.

Create a refreshing space with Refresh Labo.

* Not effective against all germs and viruses.



Refresh Labo combines three mechanisms of action for effective deodorization by transforming, containing, and preventing the generation of odorous substances.

Disinfecting Deodorizer

Refresh Labo

Clothing Fabric Air Antibacterial


<Fragrance-Free>



<Scent of the Forest>



Deodorizing Effect * The larger the enclosed area on the chart the greater the deodorizing effect.



Effectively deodorizes a variety of smells.

As well as deodorizing the air, Refresh Labo can also be used directly on floor cushions, clothing, shoe cupboards, curtains, and rugs.

Odor Elimination Study (supervised by an Olfactory Measurement Operator)

Cotton cloths with each odor applied and with and without one pump spray of each Refresh Labo were placed in individual containers, which were placed in smell-test bags. An adjudicator then smelled each smell-test bag and provided an evaluation based on the odor intensity indicator method. The chart shows the difference in odor intensity when Refresh Labo was and was not used on each odor.

The larger the difference in intensity (the larger the enclosed area on the chart) the greater the deodorizing effect. (A Niitaka study)


— <Fragrance-Free> — <Scent of the Forest>

Instructions for Use

Deodorizing an entire room	(1) Ventilate the room well. (2) Spray the entire room.
Deodorizing fabric	(1) From a distance of 20 to 30 cm, spray enough to moisten the entire area. (2) Allow to dry naturally.
Wiping down a surface/object	(1) Spray an appropriate amount on the surface/object. (2) Remove the dirt with a dish cloth.

Disinfecting Effect

Prevents the foul odors generated by bacterial growth.



Disinfection Study

Test bacterial solution was exposed to the product at room temperature for 30 seconds. After 30 seconds, the exposure was stopped, the bacterial solution was cultured in standard method agar for 48 hours at 35°C, and the bacterial count was measured. (A Niikata study)

Cleaning Effect

Walls immediately after being sprayed with Refresh Labo. Smoking residue is dripping off the wall.



Wall surface



Wall surface

! Precautions for Use

- Keep away from open flames and other sources of ignition (including static electricity).
- The appearance and labeling of this product may be altered without notice.
- Always read the "Instructions for Use" and "Precautions for Use" printed on the product before use.

Do Not Use On

Leather products, artificial leather, kimono, furs, and furniture and floors made of untreated wood or paulownia wood.

* In the case of paint, varnish, silk, rayon, other materials susceptible to water damage, and materials with waterproofing treatment, first test this product in an inconspicuous area before using over a wider area.



A dedicated spray bottle is available.

* Sold separately
Product name: Dedicated wide-neck spray bottle 500mL
Product code: 900088

Product Data

Product Name	Clothing, fabric, and air deodorizer; Store and home cleaner			
Trade Name	Refresh Labo (Fragrance-free)		Refresh Labo (Scent of the forest)	
Ingredients	Surfactant (alkylamine oxides), solvent (ethanol), disinfectant, stabilizer, metal ion sequestrants		Surfactant (alkylamine oxides), solvent (ethanol), disinfectant, stabilizer, metal ion sequestrants, fragrance	
Acidity/Alkalinity	Neutral		Neutral	
NET Volume	2.5 L	20 L	2.5 L	20 L
Packaging	2.5 L (soft bottle) × 4 units/case	20 L (Bag-in-Box container)	2.5 L (soft bottle) × 4 units/case	20 L (Bag-in-Box container)
Product Dimensions*	141 mm × 105 mm × 272 mm	294 mm × 295 mm × 287 mm	141 mm × 105 mm × 272 mm	294 mm × 295 mm × 287 mm
Case Dimensions*	293 mm × 226 mm × 289 mm	_____	293 mm × 226 mm × 289 mm	_____
Stacking Limit	3 cases	3 cases	3 cases	3 cases
Product Code	293060	293001	293160	293101
JAN Code	4975657229240	4975657229257	4975657229271	4975657229264
ITF Code	14975657229247	_____	14975657229278	_____

* Dimensions are external length × width × height. Case dimensions are score-to-score dimensions. * 20 L BiB container does not include a tap valve. Please purchase separately.